

Mattress with
mesh bedcover



3B for motion bed 90x195x8cm




聚合弹性体 床芯
Polymer Core

it can be sterilized by vapor 95~100 °C




2B for wheel chair 40x40x5cm



床芯

1H... this core be worked under 55~60 °C



1H Z-fold for mattress topper 90x188x3cm

BEING CHAMPION CO.,LTD



秉昌企業有限公司

台北市11011 信義路5段5號3樓E區27室 台北郵政12031信箱

P.O.Box: 12031, Taipei, Taiwan

Add.: Room 3E-27, No. 5, Hsin-yi Road, Sec. 5, Taipei, Taiwan 11011

Tel: +886-2-27410055 / 27016688 Fax: +886-2-27204220

e-mail: bcc@ms3.hinet.net / bcc.bcc@msa.hinet.net

Http:// bczone.net